

Surgimpex & Co.



DENTAL INSTRUMENTS CATALOGUE

History To introduce the firm briefly , M/s. Surgimpex & Co was established in 1978 Mr. Abdul Hafiz Mughal , who was well known as a specialized skilled worker and maker of Surgical Dental Veterinary Instruments.

Our Quality

- 1 Fair dealing, Straight Forwardness, Honesty, Sincerity, Truth
- 2 The Material used in the production against order is exactly as per approved samples.
- 3 Normally, in Cutting Instruments, we use (420 Type) Japanese Stainless Steel.
- 4 For non -Cutting instruments like Forceps, Tweezers etc.; we use (410 Type) Japanese Stainless Steel.
- 5 Moreover, we can also the material as per our customer requirements.

Our Price

- 1 To provide “Excellent Service and Economically Price” is the motto of our business
 - 2 All instruments are made of Stainless Steel otherwise specified.
 - 3 “All types of Mirror-Finish” “Satin/Dull-Finish” “Sand Blast-Finish” instruments are available in stock but choice of polishes will be supplied on demand without any extra charges
- Therefore, You are requested to kindly provide us a list of your desire items against which we shall appreciate to quote / offer our most competitive prices for favour of your kind approval.
Looking forward to your early kind encouraging reply.

Our Services

Manufacturing Steps

Production Capabilities

Surgimpex & Co is an employee owned company of craftsmen and has got a good reputation for quality and reliability. Our significant manufacturing steps are being described here-under briefly

1. MATERIAL USED: Always a best quality stainless steel material is purchased from steel manufacturers. Two kind of series i.e. 300, 410 and 420 of stainless steel is used for production of Surgical Instruments.

mechanization to production run

CONTENTS

01

EXTRACTING FORCEPS

Page # 03-16

02

ROOT ELEVATORS

Page # 17-32

03

**DIAGNOSTIC, PROBES, SCALERS, FILLING
& LABORATORY INSTRUMENTS**

Page # 33-83

04

**MAXILLOFACIAL SURGERY
INSTRUMENTS**

Page # 84-100

05

**ORTHODONTIC INSTRUMENTS
& PLIERS**

Page # 101-108

06


DENTAL IMPLANT INSTRUMENTS

Page#109-152

07

IMPRESSION TRAY & HOLLOW WARES

Page#153-158





SI-001
Fig. 23
Syndesmotome lower roots



SI-002
Fig. 24
Syndesmotome upper model



SI-003
Fig. 25
Syndesmotome Lower model forceps



SI-004
Fig. 79
Lower Wisdoms




SI-005
Fig. 60
Upper roots, extra slim jaws



SI-006
Fig. 62
Upper roots, extra slim jaws



SI-007
Fig. 74
Upper roots and incisors



SI-008
Fig. 33
Lower Roots



SI-009
Fig. 51
Upper Roots



SI-010
Fig. 61
Upper Roots, short extra slim jaws



SI-011
Fig. 67
Upper wisdom either sides



SI-012
Fig. 86
Lower molars decayed or with broken down crown



SI-013
Fig. 01
Upper centrals and canines



SI-014
Fig. 04
Lower incisors and canines



SI-015
Fig. 08
Lower bicusps



SI-016
Fig. 17
Upper right molars



SI-017
Fig. 73
Lower molars



SI-018
Fig. 21
Lower molars



SI-019
Fig. 02
Upper centrals and canines



SI-020
Fig. 07
Lower incisors and bicusps



SI-021
Fig. 13
Lower Canines and bicusps




SI-022
Fig. 18
Upper Left Molars



SI-023
Fig. 22
Lower molars




SI-024
Fig. 75
Lower canines and bisuspsids



Apical Set of 9 Pcs (Deep Gripe Extracting Forceps)



- No Fracture
- No Tissue Damage
- Extraction With Easy



SI-029
F-2
Lower Universal



SI-030
F-4
Lower Anterior



SI-025
F-1
Upper Universal



SI-026
F-1
Upper Universal



SI-027
F-3
Upper Anterior



SI-028
F-3
Upper Anterior

SI-031
F-4
Lower Anterior

SI-032
F-6
Upper 3RD

SI-033
F-7C
Lower Molar

A-TRAUMATIC EXTRACTION SET

Improved tapered beaks for complete root contact. Cross serrations for a never slip grasp. Easily enter the sub gingival.



Fig. 79A

Fig. 45Z



Fig. 51Z

Fig. 44Z

Fig. 46



Fig. 67

Fig. 67A

Fig. 79



SI-034

Fig. 9

Apex Forceps Lower Universal
Deep Root Pick



SI-035

Fig. 187

Apex Forceps Lower Universal
Deep Root Pick



SI-036

Fig. 189

Apex Forceps
Lower Universal
Deep Root Pick



SI-037

Fig. 190

Apex Forceps
Lower Universal
Deep Root Pick




SI-038

Fig. 46XL

Apex Forceps
Lower Root
Deep Root Pick

DEEP ROOT PICK FORCEPS (Apical Set of 9)



SI-039

Fig. 67


Lower Thrid
Molars



SI-040

Fig. 67A

Lower Thrid
Molars



SI-041

Fig. 79

Lower Thrid
Molars



SI-042

Fig. 37

Lower Root
Molars



SI-43
Fig. 2
Upper laterals and canines



SI-44
Fig. 7
Upper Premolars




SI-45
Fig. 17
Upper Molars Right



SI-46
Fig. 18
Upper Molars Left



SI-47
Fig. 67
Upper third molars



SI-048
Fig. 13
Lower Premolars



SI-049
Fig. 22
Lower Molars



SI-050
Fig. 22
Lower Molars



SI-051
Fig. 51
upper roots



SI-052
Fig. 79
Lower third molars



SI-053
Fig. 1
Upper laterals
and Canines



SI-054
Fig. 2
Upper laterals
and Canines



SI-055
73
Lower Molars



SI-056
67
Upper Third
Molars



SI-057
51A
Upper Roots



Extracting Forceps
Set of 10 Pcs
Detail of Pack
Forceps
Fig. 1
Fig. 2
Fig. 73
Fig. 67
Fig. 51A
Fig. 74
Fig. 51
Fig. 17
Fig. 7
Fig. 13



SI-058
74
Upper Roots



SI-059
51
Upper Roots




SI-060
17
Upper Roots
right



SI-061
7
Lower Premolars



SI-062
13
Lower Premolars





Set of 10 Pieces
 Fig. 1
 Fig. 7
 Fig. 13
 Fig. 17
 Fig. 18
 Fig. 22
 Fig. 33A
 Fig. 51A
 Fig. 67A
 Fig. 79



SI-063
 Fig. 1
 Upper Centrals &
 Canines



SI-064
 Fig. 7
 Upper Premolars




SI-065
 Fig. 13
 Lower Premolars



SI-066
 Fig. 17
 Upper Molars
 Right



SI-067
 Fig. 18
 Upper Molars
 Left



SI-068
 Fig. 22
 Lower Molars



SI-069
 Fig. 33A
 Lower Roots
 Narrow tips





SI-070
 Fig. 51A
 Upper Roots
 Narrow tips




SI-071
 Fig. 67A
 Upper third
 molars



SI-072
 Fig. 79
 Lower Third
 Molars





SI-080
Fig. 17 SK
Lower Molarst



SI-081
Fig. 150 SK
Upper anteriors and roots



SI-082
Fig. 151 SK
Lower anteriors and roots



SI-083
Fig. 90
Upper molars, left



SI-084
Fig. 17S
Lower, Molars



SI-085
Fig. 18LS
Upper, molars left




SI-086
Fig. 65S
Upper incisors and roots




SI-087
Fig. 69S
Upper and lower roots



SI-088
Fig. 69S
Upper and lower roots






SI-089
150
Cyber uoer incisors premolars roots



SI-090
150A
Cryer uper incisors premolars roots




SI-091
150S
Cryer upper teeth for children



SI-092
151A
Cryer upper teeth for children



SI-093
151S
Cryer upper incisors, premolars roots



SI-094
151
Cryer lower teeth for children



SI-095
286
upper incisors and roots





SI-096
89
Cohen universal forceps



SI-097
89
Upper Molars right



SI-098
23S
Lower Molars




SI-099
2
Upper and lower roots




SI-104
18R
Harris upper molars, left





SI-100
4
Upper and lower premolars




SI-105
18R
Harris upper molars, right



SI-101
7
Upper roots



SI-106
29S
Lower teeth, for children




SI-102
13
Upper incisors



SI-107
#35
Upper roots



SI-103
17
Lower molars



SI-108
Cohen
Universal forceps



SI-109
Fig. 1
Upper Anterior

SI-110
Lower Molar
Fig. 23 Cow Horn



SI-111
Fig. 74
Lower Anterior



SI-112
Fig. 88L
Upper Molar Left



SI-113
Fig. 88R
Upper Molar Right



SI-114
Fig. 74
Lower Anterior



SI-115
Fig. 106
Lower Incisors



SI-116
Fig. 174
Extra Narrow
lower Anterior
Apical Retention
Forceps



SI-117
Fig. 859
Lower Molars
Roots







SI-118 1.8ML
SI-119 2.2ML



SI-120 1.8ML
SI-121 2.2ML



SI-122 1.8ML
SI-123 2.2ML



SI-124 1.8ML
SI-125 2.2ML



SI-126 1.8ML
SI-127 2.2ML



SI-128 1.8ML
SI-129 2.2ML



SI-130
1001
Intraligamental
Anesthesia Syringe
in Black Colour



SI-131 1.8ML
SI-132 2.2ML
Angled



SI-133
Aluminium - Dark 1.8 ml



SI-134
Aluminium - Silver 1.8 ml



SI-135
Aluminium - Orange 1.8 ml



SI-136
Aluminium - Green 1.8 ml



SI-137
Aluminium - Blue 1.8 ml



SI-138
Blumenthal
15cm



SI-139
Blumenthal
15cm



SI-140
Blumenthal
15cm



SI-141
Friedman
14cm



SI-142
Friedman
14.5cm



SI-143
Micro-Friedman
14cm



SI-144
Micro-Friedman
15cm



SI-145
Micro-Friedman
16cm



SI-146
Kleinert-Kutz
18.5cm




SI-147
Beyer
18.5cm





SI-148
Mead 1A
18.5cm



SI-149
Luer
18.5cm



**HANDLES FOR ROOT ELEVATORS POINTS
A, B, C, D**



Apexo Root Elevators



Apexo



SI-150
Apexo
2.5x9mm
Straight



SI-151
Apexo
2.5x9mm
Straight



SI-152
Apexo
3.0x9mm
Straight



SI-153
Apexo
2.5x9mm
Curved



SI-154
Apexo
3.0x9mm
Curved



SI-155
Seldin
3.0x14mm
Straight

Luxator Root Elevators



Luxator



SI-156
Luxator
3mm
straight



SI-157
Luxator
3mm
curved



SI-158
Luxator
4mm
straight



SI-159
Luxator
4mm
curved



SI-160
Luxator
5mm
curved



SI-161
Luxator
5mm
curved

Bein Root Elevators



Bein



SI-162
Bein 2mm



SI-163
Bein 3mm



SI-164
Bein 4mm

Seldin Root Elevators



Seldin



SI-165
Flohr 1



SI-166
Seldin 34S



SI-167
Seldin 34

Apexo Root Elevators



Apexo



SI-168
Apexo
2.5x9mm
Straight



SI-169
Apexo
2.7x9mm
Straight



SI-170
Apexo
3.0x9mm
Straight



SI-171
Apexo
2.5x9mm
Curved



SI-172
Apexo
3.0x9mm
Curved



SI-173
Seldin
3.0x14mm
Straight

Luxator Root Elevators



Luxator



SI-174
Luxator
3mm
straight



SI-175
Luxator
3mm
curved



SI-176
Luxator
4mm
straight



SI-177
Luxator
4mm
curved



SI-178
Luxator
5mm
straight



SI-179
Luxator
5mm
curved

Bein Root Elevators



Bein



SI-180
Bein 2mm



SI-181
Bein 3mm



SI-182
Bein 4mm

Seldin Root Elevators



Seldin






SI-183
Flohr 1



SI-184
Seldin 34S




SI-185
Seldin 34



Seldin
13-3
Str.



Apical Root Elevators



SI-201
Fig. 1
str.




SI-202
Fig. 4
Left




SI-203
Fig. 5
Right


Apical Root Elevators



SI-204
Fig. 2
str.




SI-205
Fig. 3
Left



SI-206
Fig. 4
Right


Coupland Root Elevators



SI-207
Coupland No. 1




SI-208
Coupland No. 2



SI-209
Coupland No. 3

Coupland

Seldin Root Elevators



SI-210
Seldin 4L
Left



SI-211
Seldin 4L
Right



SI-212
Seldin 1L
Left



SI-213
Seldin 1R
Right




SI-214
Seldin
Left



SI-215
Seldin
Right




Lindo LeSian Root Elevators



Flohr Root Elevators



Heidbrink Root Elevators



Friedmann Roots Elevator



SI-251
Friedmann 4F
Str.

SI-252
Friedmann 5F
Left

SI-253
Friedmann 6F
Right

Heidbrink Roots Elevator



SI-254
Heidbrink
Str.

SI-255
Heidbrink
Left

SI-256
Heidbrink
Right

Wisdom Roots Elevator



SI-257
Wisdom
for 18 and 38

SI-258
Wisdom
for 28 and 48

Bayonet Roots Elevator Set



SI-259
Left


SI-260
Right

SI-261
Bayonet

SI-262
Curved.

SI-263
Straight


Luxator Roots Elevator



SI-264
Luxator
Straight

SI-265
Luxator
Curved

Roots Elevator



Luxier Roots Elevator




Luxier Roots Elevator



Claw-Lift Roots Elevator



Hylin Roots Elevator



Cryer Roots Elevator



Kopp Roots Elevator



Kopp Roots Elevator



SI-295
Kopp
Right

SI-296
Kopp
Left

Heber Roots Elevator



SI-297
Heber
Bayonet

SI-298
Heber
Bayonet


Bein Roots Elevator



SI-299
Bein
3mm

SI-300
Bein
4mm

Bein Roots Elevator



SI-301
Bein
2mm, Straight
Serrated

SI-305
Bein
2mm, Straight
Serrated

SI-309
Bein
2mm, Straight
Serrated

SI-302
2mm
Straight
Serrated

SI-303
3mm
Straight

SI-304
4mm
Curved

SI-306
2mm
Straight
Serrated

SI-307
3mm
Straight

SI-308
4mm
Curved

SI-310
2mm
Straight
Serrated

SI-311
3mm
Straight

SI-312
4mm
Curved

Seldin Roots Elevator



W. Barry Roots Elevator



Winter Roots Elevator



Pott Roots Elevator



Heidbrink Roots Elevator



SI-338 Heidbrink Straight
SI-339 Heidbrink Left
SI-340 Heidbrink Right
SI-341 Heidbrink Straight
SI-342 Heidbrink Left
SI-343 Heidbrink Right

Syndesmotomes Roots Elevator



SI-344 Straight
SI-345 Left
SI-346 Right

Periotome Roots Elevator

Syndesmotomes



SI-347 Periotome 2.1mm Straight
SI-348 Periotome 2.6mm Straight
SI-349 Periotome 2.1mm Angular
SI-350 Periotome 2mm Curved
SI-351 Periotome 2.0mm
SI-352 Periotome 2.5mm
SI-353 Periotome 3.0mm
SI-354 Periotome 4.0mm



SI-355
1S
Straight 1mm



SI-356
2S
Straight 2mm



SI-357
3S
Straight 3mm



SI-358
5S
Straight 5mm



SI-359
3C
Curved 3mm



SI-360
5C
Curved 5mm



SI-361
3CA
Curved Angel 3mm



SI-362
3IC
Inverted Curved 3mm



SI-363
DE3
Dual Edge 3mm/1.5mm



SI-364
DE5
Dual Edge 5mm/3mm

PDL ELEVATORS






SI-367
PDL Elevators
Set of 5 Pcs



Our Scalers

Compare our scalers and curettes with the best in the industry and you will immediately recognize we are equal in every way but one... price. Note the deep precision knurling for better tactile feel. The same high-grade stainless steel tips, the same light weight (under 18 grams) for less hand fatigue, the same professional craftsmanship, and surpassing the highest industry standards. Of course, any manufacturer can make these claims, but we guarantee and certify ours.



- High-polish mirror finish tips.
- Manufactured within the guidelines of the FDA quality assurance system for good manufacturing practices.
- Highest grade German stainless steel, razor sharp, durable tips. Rockwell hardened to 60 RHc.
- 27 of the most popular designs from Graceys to McCallis to Columbias, and many European patterns.
- All instruments are guaranteed against defects in material and workmanship.

- All instruments bear the CE mark, emblematic of the European community's highest standards.
- Customer satisfaction guarantee, too!
- Less than 20 grams ... Lightweight for less hand fatigue.
- Manufactured to the most exacting standards so the 50th instrument you buy is exactly like your first.
- You decide ... Why pay more for the same high-quality instruments?



TWEEZERS



SI-381
Surgical
11.5cm



SI-382
London Collage
15cm



SI-383
Semken Taylor
Str



SI-384
Semken Taylor
Cvd



SI-385
Meriam
16cm



SI-386
Dressing
14cm



SI-387
Semken Taylor
1x2 Straight



SI-388
Semken Taylor
1x2 Curved



SI-389
Gilles
15cm 1x2



SI-390
14cm 1x2



SI-391
Perry
13cm



SI-392
Atraumatic
15cm Cvd



SI-393
Atraumatic
15cm Str



SI-394
Collage with Lock
15cm Cvd



SI-395
Gerald
17.5



SI-396
With lock
15cm



SI-397
Atraumatic
17.5cm



SI-398
Adson- TC
12cm



SI-399
Adson- TC
12cm 1x2



SI-400
Adson- TC
15cm 1x2



SI-401
Adson
12cm



SI-402
Adson
12cm 1x2



SI-403
Adson
15cm 1x2



SI-404
Adson Brown
12cm



SI-405
Crane Kaplan
15cm



SI-406
Crane Kaplan
15cm



SI-407
Crane Kaplan
15cm



SI-408
Hemingway
15cm



SI-409
Miller
15cm Str



SI-410
Miller
15cm Cvd.



SI-411
Dressing Tweezer
14.5cm



SI-412
Collage Corved
Tip 15cm



SI-413
Maryam Tweezer
15cm




SI-414
Collage tweezer
15cm






SI-414-A
Perry Tweezer
12.5cm




SI-415
Self Lockingr
15cm







Pointed Probes Single Ended Handle 6mm, length 150mm



Pointed Probes Double Ended Handle 6mm, length 170mm




PERIODONTAL SET OF 10 PIECES SOLID HANDLE 6MM LENGTH 170MM



WILLIAM PROBE SET DOUBLE ENDED WITH HOLLOW HANDLE 9.5MM, LENGTH 170MM CERAMIC COATED HANDLE


Proper care must be given during cleaning, sharpening and sterilization to preserve design and integrity and to prolong usefulness. Do not expose instruments to sterilization temperatures above 350 F/177 C or to higher than recommended chemical concentrations.



PERIODONTAL INSTRUMENTS SET OF 11 PCS SOLID HANDLE 6MM, LENGTH 170MM



GRACEY SET OF 7 PCS HOLLOW HANDLE 9.5MM, LENGTH 170MM



SI-517
Gracy
1/2

SI-518
Gracy
3/4

SI-519
Gracy
5/6

SI-520
Gracy
7/8

SI-521
Gracy
9/10

SI-522
Gracy
11/12

SI-523
Gracy
13/14

GRACEY SET OF 7 PCS
HOLLOW HANDLE 9.5MM, LENGTH 170MM WITH SILICON GRIP



SI-524
Gracy
1/2



SI-525
Gracy
3/4



SI-526
Gracy
5/6



SI-527
Gracy
7/8



SI-528
Gracy
9/10



SI-529
Gracy
11/12



SI-530
Gracy
13/14

GRACEY MINI FIVE


HOLLOW HANDLE 9.5MM, LENGTH 170MM
CERAMIC COATED HANDLE



SI-531 Gracy 1/2
SI-532 Gracy 3/4
SI-533 Gracy 5/6
SI-534 Gracy 7/8
SI-535 Gracy 11/12
SI-536 Gracy 13/14

GRACEY AFTER FIVE

HOLLOW HANDLE 9.5MM, LENGTH 170MM
CERAMIC COATED HANDLE




SI-537 Gracy 1/2 AF
SI-538 Gracy 3/4 AF
SI-539 Gracy 5/6 AF
SI-540 Gracy 7/8 AF
SI-541 Gracy 11/12 AF
SI-542 Gracy 13/14 AF

Proper care must be given during cleaning, sharpening and sterilization to preserve design and integrity and to prolong usefulness. Do not expose instruments to sterilization temperatures above 350 F/177 C or to higher than recommended chemical concentrations.

HOLLOW HANDLE 10MM, LENGTH 170MM

Curetters Set of 8 Pieces



SI-543
Curette
Columbia
2L-2R

SI-544
Curette
Columbia
4L-4R

SI-545
Curette
McCall
13S-14S

SI-546
H6-H7
Scaler
Sickle


SI-547
H5-U15
Scaler
Sickle

SI-548
H5-13
Scaler
Sickle

SI-549
H15-30
Scaler
Sickle

SI-550
H15-33
Scaler
Sickle

FILLING INSTRUMENTS (Gold Pointed)
SOLID HANDLE 6MM, LENGTH 170MM




FILLING INSTRUMENTS
spatula 6mm solid handle
length 170


FILLING INSTRUMENTS
6mm
length 175mm

FILLING INSTRUMENTS
ball 6mm
length 170mm

FILLING INSTRUMENTS
ball spatula
6 mm
length 175mm


FILLING INSTRUMENTS
plan 6 mm
length 165mm






ENDODONTICS INSTRUMENTS


ANTERIOR



LUKS








COMPOSITE/PLASTIC FILLING INSTRUMENTS SOLID HANDLE 6MM, LENGTH 170MM

BLACK


PLUGGERS/CONDENSERS

Hollenback





PLUGGERS/CONDENSERS

RESTORATIVE INSTRUMENTS SOLID HANDLE 6MM, LENGTH 170MM




PLUGGERS/CONDENSERS





RESTORATIVE INSTRUMENTS SOLID HANDLE 6MM, LENGTH 170MM

BURNISHERS




PERIODONTAL FILES




COMPOSITE/PLASTIC FILLING INSTRUMENTS SOLID HANDLE 6MM, LENGTH 170MM

BURNISHER FILLING INSTRUMENTS




UNIVERSAL FILLING INSTRUMENTS




FILLING INSTRUMENTS SOLID HANDLE 6MM, LENGTH 170MM






AMALGAM INSTRUMENTS SOLID HANDLE 6MM, LENGTH 170MM

Excavators Instruments Set of 11 Pieces



SICKLE SCALERS INSTRUMENTS SOLID HANDLE 6MM, LENGTH 170MM


SICKLE SCALERS




PERIODONTAL INSTRUMENTS SOLID HANDLE 6MM, LENGTH 170MM


JACQUETTE SCALERS

CRANE-KAPLAN




PERIODONTAL INSTRUMENTS SOLID HANDLE 6MM, LENGTH 170MM SCALER/CURETTE SETS





GRACEY CURETTE INSTRUMENTS SOLID HANDLE 6MM, LENGTH 170MM GRACEY CURETTES



SI-757
9/10t
(Shown)
#2,#4,#6,
#7,#8,#9

SI-758
9/10
Rigid
#2,#4,#6,
#7,#8,#9

SI-759
11/12
Rigid
#2,#4,#6,
#7,#8,#9

SI-760
11/12t
(Shown)
#2,#4,#6,
#7,#8,#9

SI-761
11/12
Extra Rigid
#4,#6,#7,
#8,

SI-762
13/14t
(Shown)
#2,#4,#6,
#7,#8,#9

SI-764
13/14
Extra Rigid
#4,#6,#7,
#8,

SI-765
15/16t
(Shown)
#2,#4,#6,
#7,#8,#9

SI-766
15/16
Rigid
#2,#4,#6,
#7,#8,#9

SI-767
17/18
#2,#4,#6,
#7,#8,#9

SI-768
17/18
Rigid
#2,#4,#6,
#7,#8,#9

CURETTE INSTRUMENTS
HOLLOW HANDLE 8.5MM, LENGTH 170MM



SI-769
Hemingway
2mm,2mm

SI-770
Hemingway
3mm,3mm

SI-771
Hemingway
4mm,4mm

CURETTE INSTRUMENTS
SOLID HANDLE 6MM, LENGTH 170MM



SI-772
Hemingway
2mm,2mm

SI-773
Hemingway
3mm,3mm

SI-774
Hemingway
4mm,4mm

CURETTE INSTRUMENTS
SOLID HANDLE 6MM, LENGTH 170MM





SI-775
Hemingway
2mm,2mm

SI-776
Hemingway
3mm,3mm

SI-777
Hemingway
4mm,4mm










EXCAVATORS INSTRUMENTS SOLID HANDLE 6MM, LENGTH 170MM

ENGLISH PATTERNS









HOES


CARVERS

AMALGAM FILES




RESTORATIVE INSTRUMENTS SOLID HANDLE 6MM, LENGTH 170MM


MARGIN TRIMMERS



P.K. THOMS WAXING INSTRUMENTS SOLID HANDLE 6MM, LENGTH 170MM

P.K. THOMAS WAXING INSTRUMENTS







SI-893
00

SI-894
0

SI-895
1

SI-896
2

SI-897
3

SI-898
4

SI-899
3
Weston

SI-900
1
medium cement

SI-901
1



SI-902
Le Corn

SI-903
Le Corn

SI-904
Le Corn

SI-905
MSY

SI-906
MSY

SI-907
Le Corn

SI-908
Hyllin

SI-909
Hyllin

SI-910
Hyllin

SI-911
Zahle

SI-912
Zahle

SI-913
Zahle

SI-914
Gritman

ALGINATE, PLASTER SPATULAS, WAX KNIVES PLASTER FOR KNIVES & PLASTER SHEAR



SI-915
Wooden Handle
Fig. 5
21.5cm



SI-916
Wooden Handle
Fig. 7
19cm



SI-917
Wooden Handle
Fig. 8
20cm curved



SI-918
Wooden Handle
Fig. 9
222.5cm curved



SI-919
Hook Knives



SI-920
Fig. 1
Lessman
12.5cm



SI-921
Fig. 2
Lessman
17cm



SI-922
Fig. 1
with large spoon
12.5cm



SI-923
Fig. 2
with large spoon
17cm



SI-924
Bacolite Handle
Wax Knives
Small



SI-925
Bacolite Handle
Wax Knives
Large



SI-926
Wooden Handle
Wax Knives
Small



SI-927
Wooden Handle
Wax Knives
Large



SI-928
Plaster Knife



SI-929
Fig. 1
17 cm



SI-930
Fig. 3
16 cm



SI-931
Fig. 4
16 cm



SI-932
Special Inox
19.5 cm

**SI-933**

Technical Datas:

Intercondylar distance	120 mm.
condylar inclination	fixed at 30°
Bennet angle	fixed at 7°
Incisal table inclination	10°
Distance between mounting plates	85 mm.
Weight	400 gr.
Opening (item 5000)	125°
Opening (item 5010)	180°

**SI-934**

Technical Datas:

Intercondylar distance	110 mm.
condylar inclination	fixed at 30°
Bennet angle	fixed at 7°
Incisal table inclination	flat
Distance between mounting plates	85 mm.
Weight	450 gr.
Opening	180°

**SI-935**

Stainless Steel Finishing

SI-936
Royal Blue

**SI-937****SI-938****SI-939****SI-940****SI-941****SI-942****SI-943****SI-944****SI-945**

SI-946
Nakawalni
Dental Laboratory
Instruments



SI-947
Mouth Gag Molt
11cm for Children



SI-948
Mouth Gag Molt
14cm for Adult



Cheek Retractor
SI-949 140 m
SI-950 180 m



Cheek Retractor
SI-951 140 m



SI-952
Vestibule
15 cm with or without
lip protection



SI-953
Wangen Holder
15 cm with or without
lip protection



SI-954
Wangen Holder
with dual function 17cm



SI-955
140 mm

SI-956
150 mm

Cheek Retractor



SI-957
Fig. 1
Weitlaner
10.5cm
13.5cm



SI-958
Fig. 2
Weitlaner
10.5cm
13.5cm



SI-959
Fig. 2
Stainless Steel
Double
Ended Cheek
REtractors



SI-960
Fig. 2
Stainless Steel
Single
Ended Cheek
REtractors



Weitlander
SI-961 11 cm
SI-962 13 cm
SI-963 15 cm



Jennings
SI-964 11 cm
SI-965 13 cm
SI-966 15 cm




SI-967
Denhard
12 cm



SI-968
Rose Konig
16 cm
12 cm



SI-969
Rose Konig
19 cm





SI-970
Cheek Retractor
Long



SI-971
Cheek Retractor
Long



SI-972
Cheek Retractor
Long



Retractor with
Ergo Grip
SI-973 24cm
ergonomic handle,



SI-974
Cheek Retractors,
Retractor Ivory
with adjustable
& fixable, 20 cm



Cheek Retractors,
tongue depressors
SI-975 19cm



Cheek Retractors,
Minnesota
SI-976 16cm



Cheek Retractors,
SI-977 20cm
easy to bend
by hand



Cheek Retractors,
retractor and tongue
pushers BISHOP with
comfortable margin
SI-978 16cm



Cheek Retractors,
retractor combination
combination
SI-979 25cm



SI-980
Retractors, Unique
Length 175mm
angled 45



SI-981
Retractors, Unique
Length 175mm
angled 90



SI-982
Retractors, Unique
Length 175mm
angled 45



SI-983
Retractors, Unique
Length 175mm
angled 90



SI-984
Retractors
16cm



SI-985
Retractors
19cm



SI-986
Retractors, Mini-Langenbeck
16cm, hooks 22mm



SI-987
Retractors, Cheek
Retractor Combination
21cm, 2



SI-988
Mini Retractors
about 8cm



SI-989
Retractors
17cm, hooks 22mm



SI-990

Retractors, surgical-flat
surgically-flat, 21cm, hook
30x16 mm



SI-991

Retractors, surgical-flat
surgically-flat, 21cm, hook
42x11 mm



SI-992

Retractors, surgical-flat
surgically-flat, 21cm, hook
63x12 mm



Retractors, surgical-flat
surgically-flat, 21 cm, hook

SI-993 30x16mm
SI-994 42x11mm
SI-995 63x12mm



Retractors, Langenbeck
21cm, hook

SI-996 30x11 mm
SI-997 30x16 mm



SI-998

Retractors
21 cm



SI-999

Farabeuf
(120 x 12mm
Set of 2 pices



SI-1000

Retractor Seldin



SI-1001

FELDMAN
RETRACTORS
5-1/2", 14cm



SI-1002
Iris
115mm
Straight



SI-1003
Iris
115mm
Curved



SI-1004
Goldman-Fox
130mm
Saw edge, straight



SI-1005
Goldman-Fox
130mm
Saw edge, curved



SI-1006
La Grange
115mm
Saw edge, double curve



Mayo
SI-1007 140mm
SI-1008 170mm
SI-1009 230mm
Straight



Mayo
SI-1010 140mm
SI-1011 170mm
SI-1012 230mm
Curved



Metzenbaum
SI-1013 140mm
SI-1014 180mm
SI-1015 200mm
Straight



Metzenbaum
SI-1016 140mm
SI-1017 180mm
SI-1018 200mm
Curved



SI-1019
Kelly
160mm
Straight



SI-1020
Kelly
160mm
Curved



SI-1021
Iris
115mm
Straight



SI-1022
Iris
115mm
Curved



SI-1023
La Grange
115mm
Double Curved



SI-1024
Goldman-Fox
130mm
Straight



SI-1025
Goldman-Fox
130mm
Curved



SI-1026
Kelly
160mm
Curved



SI-1027
Iris
115mm
Straight



SI-1028
Iris
115mm
Curved



SI-1029
Micro Iris
90mm
Curved



SI-1030
Iris
115mm
Angled

SI-1031
Iris
115mm
Curved

SI-1032
Goldman-Fox
130mm
Saw edge,
straight

SI-1033
Goldman-Fox
130mm
Saw edge,
curved

SI-1034
Goldman-Fox
130mm
Saw edge,
curved

SI-1035
La Grange
115mm
Saw edge,
double
curved

SI-1036
Neumann
120mm
Saw edge,
double
curved



SI-1037
Kilner
120mm
Straight

SI-1038
Kilner
120mm
Curved

SI-1039
Joseph
140mm
Straight

SI-1040
Joseph
140mm
Curved

SI-1041
Kelly
160mm
Straight

SI-1042
Kelly
160mm
Straight

SI-1043
Locklin
160mm 'Saw edge
angular

SI-1044
Locklin
120mm
160mm 'Saw edge
angular

SI-1045
Dean
170mm 'Saw edge
angular



SI-1046
Spencer
90mm

SI-1047
Spencer
130mm

SI-1048
Spencer
115mm

SI-1049
Noyes
110mm, straight

SI-1150
Noyes
110mm,
curved

SI-1151
Noyes
110mm,
Angular

SI-1152
CastroSlej
150mm,
curved



SI-1053
130mm
Straight



Mayo
SI-1054 130mm
SI-1055 160mm
SI-1056 280mm
Straigh Curved



Mayo
SI-1057 140mm
SI-1058 200mm
Straigh Curved



Metzenbaum
SI-1059 145mm
SI-1060 200mm
Straigh Curved



SI-1061
Beebee
110mm
Saw edge
straight




SI-1062
Beebee
110mm
Saw edge
curved



SI-1063
Beebee
110mm
Saw edge
straight



SI-1064
Beebee
110mm
Saw edge
curved



SI-1065
Universal
120mm
Saw edge



SI-1066
Universal
120mm



Crile
SI-1067 130mm
SI-1068 150mm
SI-1069 180mm



SI-1070
Halsey
130mm



SI-1071
Baumgartner
130mm



SI-1072
Derf
120



SI-1073
Delicate
150mm



Mayo-Hegar
SI-1074 130mm
SI-1075 160mm
SI-1076 180mm
SI-1077 200mm



SI-1078
Crile-wood
150mm
straight



SI-1079
Crile-wood
150mm
Curved



Castroveijo
SI-1080 130mm
SI-1081 150mm
SI-1082 180mm



SI-1083
 Halsey
 130mm



SI-1084
 Baumgartner
 130mm



SI-1085
 Derf
 120mm



SI-1086
 Delicate
 150mm



SI-1087
 Adson
 170mm



Mayo-Hegar
SI-1088 130mm
SI-1089 150mm
SI-1090 180mm
SI-1091 200mm



SI-1092
 Crile wood
 150mm



Castroveijo
SI-1093 str 140mm cvd
SI-1094 str 180mm cvd



Olsen-Hegar Needle
 Holder with scissors
SI-1095 140mm



Gillies Needle
 Holder with scissors
SI-1096 160mm



Mathieu
SI-1097 140mm
SI-1098 180mm



Castroveijo
SI-1099 str 140mm cvd
SI-1100 str 180mm cvd



Olsen-Hegar Needle
 Holder with scissors
SI-1101 140mm



Gillies Needle
 Holder with scissors
SI-1102 160mm




Hartmann Baby Mosquito
SI-1103 90mm straight



SI-1104
 125mm
 Curved

SI-1106
 Curved

Halstrad Mosquito
SI-1105 125mm straight





SI-1143
Rubber Dam



SI-1144
Anisworth



SI-1145
Brewer



SI-1146
Ivory light weight



SI-1147
Palmer



SI-1148
Seaparing Elastic
Plier



SI-1149
Fig. 15



SI-1150
Fig. 18



SI-1151
Fig. 24



SI-1152
Fig. 25



SI-1153
Fig. 26



SI-1154
Fig. 28



SI-1155
Fig. 30



SI-1156
Fig. 31



SI-1157
Fig. 51



SI-1158
Fig. 138



SI-1159
Fig. 139



SI-1160
Fig. 200



SI-1161
Fig. 201



SI-1162
Fig. 202



SI-1163
Fig. 203



SI-1164
Fig. 204



SI-1165
Fig. 205



SI-1166
Fig. 3



SI-1167
Fig. W3



SI-1168
Fig. W3



SI-1169
Fig. 56A



SI-1170
Fig. 7



SI-1171
Fig. 7A



SI-1172
Fig. 8



SI-1173
Fig. 8A



SI-1174
Fig. W8A



SI-1175
Fig. 10



SI-1176
Fig. 11



SI-1177
Fig. 12A



SI-1178
Fig. 13A



SI-1179
Fig. 14



SI-1180
Fig. 14A



SI-1181
Fig. 14AN



SI-1182
Fig. 22



SI-1183
Fig. 27



SI-1184
Fig. 207



SI-1185
Fig. 29



SI-1186
Fig. 206



SI-1187
Fig. 208



SI-1188
Fig. 209



SI-1189
Fig. 1



SI-1190
Fig. 1A



SI-1191
Fig. 2



SI-1192
Fig. 2A



SI-1193
Fig. W2A



SI-1194
Fig. 0



SI-1195
Fig. 00



SI-1196
Fig. 210



SI-1197
Fig. 211



SI-1198
Fig. 212



SI-1199
Fig. 6



SI-1200
Fig. 9



SI-1201
Fig. N06



SI-1202
Fig. U67



SI-1203
Fig. 214



SI-1204
Fig. P1



SI-1205
Fig. P2



SI-1206
Fig. G-1



SI-1207
Fig. G-2



SI-1208
Power Crown Remover Lift



SI-1209
Power Crown Remover Right



SI-1210
Crown Speraders



SI-1211
Mershon
14cm
Band Pusher




SI-1212
Fig. 0



SI-1213
Manual Crown remover



SI-1214
Complete with 3 tips
and tightening keys



Semi automatic crown



SI-1215
Manual Crown Remover



SI-1217
Fig. 5
Ivory 8
5 cm



SI-1218
Fig. 6
Ivory 8
6 cm



SI-1219
Fig. 7
Ivory 8
7 cm



AE-1220
Fig. 1
Tofflemire Universal

Sigveland Retainers



AE-1221
Fig. 2
Tofflemire Junior



AE-1222
Fig. 3
Tofflemire Senior



SI-1223
Fig. 3
Narrow



SI-1224
Wide

Amalgam Well



SI-1225
Stainless Steel
Amalgam Well



SI-1226
Ivory 12



SI-1227
Ivory 1



SI-1228
Ivory



SI-1229
Elliot



SI-1230
Fig. 10
165
English Pattern



SI-1231
Fig. 1
185mm
German Pattern



SI-1232
Fig. 1
150mm
Swiss Pattern



SI-1233
Fig. 1
150mm



SI-1234
Fig. 2
150mm, 45



SI-1235
Fig. 3
140mm,



SI-1236
Fig. 3
5" (12.7cm)



SI-1237
Fig. 3G
5" (12.7cm)



SI-1238
Fig. 4
5" (13.7cm)



AE-1239
6 way Adjustable Angled
Scalpel Handle

AMALGAM CARRIERS





SI-1252
Karaner NeSIns
Bone Chisel



SI-1253
Karaner NeSIns
Bone Chisel



AE-1254
Molt 18.5cm
Auto Clavable
Plastic Tips



AE-1255
Molt 18.5cm
Non Auto Clavable
Plastic Tips



SI-1256
Miller



SI-1257
Miler

BONE FILES



SI-1258
Miler 21



SI-1259
Miler 1X
cross cut



SI-1260
Miler 2X
cross cut



SI-1261
Miler 52



SI-1262
Wahl 1

PERIOSTEEL ELEVATORS



SI-1263
BUSER



SI-1264
Ohl 10



SI-1265
Howard 9H



SI-1266
Freer



SI-1267
Freer 15



SI-1268
Freer 16



SI-1269
Molt 9



SI-1270
Allen 9A



SI-1272
Seldin 23



SI-1273
Goldman fox
14



SI-1274
P24G



10.0mm

SI-1275
P24G



6.4mm

SI-1275A
Raspator Raspatory
Length 175mm



3.0mm

SI-1275B



4.0mm

SI-1276



5.0mm

AE-1276A
Raspator Willig
Length 170 mm Differ



Behind Every Smile is a Story

QUALITY - ЯКІСТЬ • PRECISION - ТОЧНІСТЬ • ACCURACY - ПРАВИЛЬНІСТЬ • DURABILITY - ЗНОСОСТІЙКІСТЬ




STYLE A

STYLE B

- Super-sharp diamond bevel cutters
- smooth ultra-polished and contoured edges
- brilliant cross-hatched stainless steel finish
- superior T-12 stainless steel end inserts
- sliding lock joint for balance & precise alignment



SI-1277
Aluminium Height
Gauge 0.22



SI-1278
Aluminium Height
Gauge 0.18



SI-1279
Bone Bracket Gauge



SI-1280
Bone Thickness Gauge
40mm 180mm



SI-1281
Measuring Forceps for Jaw
Ridge width measurement



SI-1282
Direct bonding
bracket holders
135mm



AE-1283
Delux Bracket
Holder Tweezers



SI-1284
Direct bonding
bracket holders
136mm



SI-1285
Direct bond
bracket holders
tweezer
Small tip



SI-1286
Direct bonding
bracket holders
140mm



SI-1287
Direct bonding
bracket holders
135mm



SI-1288
Direct bonding
bracket holders
140mm



SI-1289
Band seater auto cleavable



SI-1290
Ortho light scaler double ended



SI-1291
Ortho light scaler double ended



SI-1292
Ortho light scaler double ended



AE-1293
Twiel on elastic module instruments



AE-1294
Twiel on elastic module instruments



SI-1295
Schure Band Seater



SI-1296
Ortho light Seater



SI-1297
Ortho light Seater



SI-1298
Ortho light Seater heavy



SI-1299
Ortho instruments



SI-1300
Ortho instruments



SI-1301
Ortho Light Band Pusher Scaler



SI-1302
Ortho light scaler double ended
ligature director



SI-1303
Ortho light scaler double ended
ligature director



SI-1304
Ortho light scaler double ended
ligature director scaler director tivate coil spring



SI-1305
Ortho band pusher scaler



SI-1306
Ortho light band pusher scaler



SI-1307
Ortho lighture director single ended



SI-1308
Ortho light scaler double ended
ligature director scaler directod



SI-1309
Ortho light scaler double ended
ligature director scaler director tivate coil spring



SI-1310
Ortho light double ended
ligature director



SI-1311
Ortho lighture director scaler



SI-1312
Ortho light lighture director double ended



SI-1313
double ended cement spatula



SI-1314
Ortho light scaler single ended


BOX JOINT PLIERS

1315
Universal
Wire Bender



SI-1315
Universal Dolphin Style
wire bending & cutting pliers

SI-1321
Bracket RemoSing Plier curved
eyes removal of brackets in the
posterior region Length 130mm



SI-1316
The curved blade of the band
remoSing plier can strip easily
and gently remove to
the teeth length 130 mm



SI-1322
young wire bending &
Loop forming pliers



MC-1317
Tweed Loop
forming pliers



SI-1323
wire bending &
Loop forming pliers



SI-1318
curved plier
Length 130mm



SI-1324
Adams small wire bending
and loop forming plier



SI-1319
curved plier
Length 135mm



SI-1325
Bracket remoSing Plier
Length 135mm



SI-1320
Nance Loop
forming Pliers



SI-1326
Weingart Plier
Length 145mm



BOX JOINT PLIERS

SI-1327
Adams Plier
Length 125mm



SI-1333
Adams Plier
Length 125mm



SI-1328
Falcon Beak
contouring pliers



AE-1334
Tweed
loop forming pliers



SI-1329
wires up to 0.020" 0.7mm



SI-1335
Bird Beak Plier Standard
Length 125mm



SI-1330
Johnson Coutouring
Pliers



SI-1336
Mini Ligature cutters



SI-1331
Large wire bending and
loop forming pliers



SI-1337
Small
Fluch cut distal and cutting pliers



SI-1332
Torquing Pliers



SI-1338
Large
Fluch cut distal and cutting pliers



BOX JOINT PLIERS

AE-1339
Tweed Angel
1.78mm



SI-1345
Contouring Plier



AE-1340
Tweed Angel
1.78mm



SI-1346
Bands & crowns plier
in order to skirt




SI-1341
Lingula/ Palatal arch
forming pliers



SI-1347
Ligature forming Pliers



SI-1342
Aderer
facebow loop forming pliers



SI-1348
length 0.91mm



AE-1343
Waldsavs wire bending
and cutting pliers



SI-1349
Bird beak wire bending
& cutting Pliers



SI-1344
Aderer Plier
Length 130mm



SI-1350
Weingart Plier medium
Length 145mm



LAP JOINT PLIERS



SI-1351
Distal End Cutter
Regular With T.C



SI-1352
Hard Wire Cutter With T.C




SI-1353
Ligature Cutter Angled
With T.C



SI-1354
Light Wire Cutter Micro
With T.C




SI-1355
Reversed Ligature Cutter Angled
With T.C



SI-1356
Posterior Band Remover Pliers
Long With T.C




SI-1357
Posterior Band Remover Pliers
Medium With T.C




SI-1358
bracket Remover Pliers
Straight With T.C



SI-1359
Bird Break Plier With T.C



SI-1360
Light Wire Bird Beak
Long Plier With T.C




SI-1361
Adam Plier With T.C



SI-1362
Tweed Arch Forming
Plier With 1.78mm T.C



SI-1363
Weingart
Plier Slim With T.C




SI-1364
Bracket Remover
Plier Angled With T.C

LAP JOINT PLIERS



SI-1365
Weingart
Plier Regular With T.C



SI-1372
Posterior Bend
Remover Plier
Long




SI-1366
Weingart General Utility Remover
Plier Long Mouth With T.C




SI-1373
Posterior Bend Remover
Plier Short



SI-1367
How Plier Straight With T.C



SI-1374
Posterior Bend Remover
Plier Medium



SI-1368
How Plier Curved With T.C




SI-1375
Bracket Remover
Plier Straight




SI-1369
Step Plier Angled With
T.C 0.55mm



SI-1376
Bracket Remover
Plier Angulated



SI-1370
Weingart Plier Utility With T.C



SI-1377
How Plier Straight



SI-1371
Weingart Plier Universal 45
With T.C



SI-1378
Bird Beak Plier
Angled Bending

LAP JOINT PLIERS



SI-1379
How Plier Curved



SI-1386
Tweed Cutter Plier
With Groove



SI-1380
Hammer Head Tie Back Plier




SI-1387
Round Nose Plier



SI-1381
Arch Forming Plier



SI-1388
Adam Plier



SI-1382
Regular Tweed Arch Forming
Plier 1.78mm



SI-1389
Jara Break Plier




SI-1383
Bird Beak With Cutter Pliers
3 Groove



SI-1390
Band Seating Plier



SI-1384
Light Wire Bird Beak Plier
With Groove



SI-1391
Light Wire Bird Beak
Long Plier



SI-1385
Tweed Loop Plier



SI-1392
Light Wire Plier Two
Step Groove
0.51mm/.020"

LAP JOINT PLIERS




SI-1393
Ball Hook Crimping Plier



SI-1394
Marcotte Looping Plier



SI-1395
Tweed Loop Plier



SI-1396
Palatal Arch Plier



SI-1397
Nance Plier



SI-1398
Younger Plier With Grooved



SI-1399
Wingurt Plier Slim




SI-1400
Tweed Loop Forming Plier TC



SI-1401
Wingurt Utility Plier



SI-1402
De La Rosa Plier
With Grooved



SI-1403
De La Rosa Plier



SI-1404
Crown And Band
Contouring Plier



SI-1405
Flat Nose Plier With
Serration



SI-1406
Flat Nose Plier Smooth



SI-1407
Stop V Bending Plier

LAP JOINT PLIERS



SI-1408
Hollowchop Contouring Plier



SI-1409
Lap-Joint Three Prong Plier




SI-1410
Three Prong Pliers Regular



SI-1411
Weingart Plier



SI-1412
Lingual Plier Long Handle
90 Angled



SI-1413
Lingual Weingart
Fine Long Handle
60 Angled



SI-1414
Direct Bonding Bracket
Remover Plier TC



SI-1415
Pin & Ligature Cutter TC



SI-1416
Mini Pin & Ligature Cutter TC



SI-1417
Standard Ligature Cutter TC



SI-1418
V Bend Stop Plier TC



SI-1419
Angled Direct Bond Bracker
Remover TC

SI-1420

Distal end cutter safety hold plier



SI-1421

Hook crimping plier Angled
Straight TC



SI-1422

3 prong clasp adjustable
plier TC



SI-1423

Utility Plier Weingart TC



SI-1424

Utility Plier Weingart TC



SI-1425

Straight Direct Bond Bracket
Remover TC.



SI-1426

Pin Holding Plier TC.



SI-1427

Hard wire cutters



SI-1428

Elastic Separating Plier
Round Handle



SI-1429

Elastic Separating Plier



SI-1430

SI-1431

Smaha
Elastomeric
ligature
placing forceps
14mm



SI-1432

Elastomeric
ligature
placing forceps
with groove
falcon 140mm




SI-1433

Elastomeric
ligature
placing forceps
with groove
falcon 140mm



Normal Anatomy – Dental Implants

Bone Level Normal
(No Bone Loss)



IMPLANT GAGES



IMPLANT GAGES





SI-1444
Clamps Osteometer
Length 155 mm
range 0-25mm



SI-1445
Measuring Loop
Straight length 90mm
measuring range
0-40mm



SI-1446
Measuring Loop
angled length 90mm
measuring range
0-40mm



SI-1447
Measuring Loop
CastroSlejo
angled length 175 mm
measuring range
0-40mm



SI-1448
Boley Gauge
Stainless steel



SI-1449
Measuring Loop
angled length
175 mm
measuring range
0-40 mm



SI-1450
Measuring Loop
Straight length
175 mm
measuring range
0-40 mm



SI-1451
Fig. 1
Iwanson
spring caliper
for metal and
porcelain



SI-1452
Fig. 2
Iwanson
spring caliper
for wax




SI-1453
Bone thickness Gauge
40mm, 180mm



SI-1454
Gauge
for measuring
distance
0-90mm

BLACK TITANIUM COLOUR COATED INSTRUMENTS



SI-1455
#1

SI-1456
#2

AE-1457
#3
extra flex

SI-1458
Micro mini

SI-1459
Mini 1

AE-1460
Mini 4
Extra flex


AE-1461
4
Extra flex

AE-1462
#5
flex thin

AE-1463
#6
flexi thin

AE-1464
Mini 3
Extra flex

BLACK TITANIUM COLOUR COATED INSTRUMENTS



SI-1465
Ab1
Boghosian

SI-1466
Ab2
Boghosian

SI-1467
W3

SI-1468
8A

SI-1469
4/5
Gregg

SI-1470
Interproximal
Carver


SI-1471
Interproximal
Carver Long

SI-1472
A6

SI-1473
Contouring
instruments

BURNISHERS

FREEDMAN BURNISHERS



SI-1474
21B

SI-1475
27/29

SI-1476
2
Ladmore

SI-1477
3
Ladmore


SI-1478
Freedman
Duckhead

SI-1479
Freedman
Small Contact
Forming

SI-1480
Freedman
Large Contact
Forming

SI-1481
Small
Medium Ball
Burnisher

BLACK TITANIUM COLOUR COATED INSTRUMENTS



SI-1482
A
Cosmetic
Contouring

SI-1483
B
Cosmetic
Contouring

SI-1484
C
Cosmetic
Contouring

SI-1485
D
Cosmetic
Contouring

SI-1486
E
Cosmetic
Contouring


SI-1487
F
Cosmetic
Contouring

SI-1488
Marginal Ridge
& Embrasure
Shaping
Instruments


SI-1489
Occlusal
Anatomy
Instruments

SI-1490
Composite
Packing
Instruments

BLACK TITANIUM COLOUR COATED INSTRUMENTS



SINUS LIFT INSTRUMENTS HOLLOW HANDLE 9.5MM LENGTH 170MM



SINUS LIFT INSTRUMENTS HOLLOW HANDLE 9.5MM LENGTH 170MM



SI-1506
Sinus Lift Elevator



SI-1507
Sinus Lift Elevator



SI-1508
Sinus Lift Elevator



SI-1509
Sinus Lift Elevator



Sinus Lift Instrument,
Flexible Sinus Lift instruments
reduce the risk of perforation of the
"Schneider's membrane"



SI-1526
Fig. 1

SI-1527
Fig. 2

SI-1528
Fig. 3

Sinus Lift Elevator,
10 mm Straight grip
Length 175 mm, Sinus lift elevators, dull,
for gentle peeling and lifting the
Schneider membrane



SI-1529
Fig. 1

SI-1530
Fig. 2

SI-1531
Fig. 3

Sinus lift instrument
10 mm Straight grip Length 180 mm,
Buser (B1, B2, B3) raspatories,
quality medical stainless steel,
corrosion-resistant




SI-1532
B01

SI-1533
B02


SI-1534
B03

Sinus lift instrument 10 mm Straight grip Sinus lift elevators, blunt or sharply curved
(H01, H02, K03, K04, K78, K77)



SINUS LIFT INSTRUMENT

Sinus Lift Elevators, 10 mm Hollow Handle length 170mm Blunt Curved Set of 6 Pcs



SI-1541

SI-1542

SI-1543

SI-1544

SI-1545

SI-1546

IMPLANT SINUS LIFT INSTRUMENTS
PERIODONTAL SURGERY SURGICAL CHISELS
KRAMER NEVINS SINUS LIFT



SI-1547 TG



SI-1548 Solt3/4



SI-1549 TGO



SI-1550 29K



SI-1551 FEDI4



SI-1552 OCH3



SI-1553 OCH4



SI-1554 FEDI1



SI-1555 FEDI2



SI-1556 FEDI3



SI-1557 Gf36

IMPLANT

SURGICAL CHISELS



SI-1558 13K/TG



SI-1559 13K/TG Titanium



SI-1560 13K/TG



AE-1561 13K/TG Titanium



AE-1562 13K/KL 4/6mm
Titanium Available



AE-1563 Pt2 Titanium
A/V Regular



AE-1564 Pt1 Titanium
A/V Regular



AE-1565 Pt3 Titanium
Angel S/E



AE-1566 Pt4 Titanium
STR S/E



SI-1567 Pt5 Titanium
TG CVD S/E

MIXING SPATULA WITH SHOVEL



AE-1568
Bone Carrier
Titanium



SI-1569
Raspator Prichard
Length 180 mm
10 mm hollow handle



SI-1570
Raspator Prichard
Length 180 mm
10 mm hollow handle
double ended



SI-1571
Raspator Prichard
Length 180 mm
10 mm hollow handle
double ended
10.5mm 4.5mm



SI-1572
Mixing Spatula
with shovel
17 cm



SI-1573
Mixing Spatula
with shovel
18.5 cm

BONE GRAFTING INSTRUMENTS COMPACTOR /TAMPERS

Product offering includes instruments for harvesting, mixing and compacting bone Provides clinicians with the proper instrumentation for bone grafting procedures Weighted bone basin Prevents tipping and spilling of bone and other materials during mixing



AE-1574
Graft Packer
Concave Tips
DE 4/6mm

SI-1575
Bone Spoon
Packer DE 3mm

AE-1576
Graft Packer
Serrated Tips
DE 6/8mm

AE-1577
Graft Packer
Serrated Tips
DE 3/4mm

AE-1578
Graft Packer
Smooth Tips
DE 5/8mm

SI-1579
titanium Graft Packer
placement site depth
gauge 3.5 mm

SI-1580
Graft Material
Convave tip for large Gr
placement site depth
gauge 3.5 mm

BONE GRAFTING INSTRUMENTS



SI-1581
Blade and Pusher
Length 180 mm
for straight or angled
blade 15x5 mm
Pushser



SI-1582
Blade and Pusher
Length 180 mm
for straight or angled
blade 15x5mm
Pusher



SI-1583
Spoon and Pusher
Length 180 mm
for straight or angled
blade 10x5mm
Pusher



SI-1584
tunnelierungs
instruments
length 175 mm



MC-1585
tunnelierungs
instruments
length 175 mm
different designs



SI-1586
Tammer
Length 170mm
for bone material
double ended



BONE GRAFTING INSTRUMENTS

1587 4mm
1588 5mm
Osmactor
150mm

SI-1589
Sinus Window
Bone Tampers
4mm

SI-1590
Sinus Window
Bone Tampers
7mm


SI-1591
6/8mm

SI-1592
3/4mm

Pluggers Clarizio
Length 160 mm, with titanium tips for mixing and pouring of autologous, heterologous bone substitute material, diameter (3/4 or 6/8 mm)

SPOON AND PUSHER (TRINOVO)

Trinovo handle, different workers, for bone material, high-grade medical stainless steel, corrosion-resistant



SI-1593
10x5mm

SI-1594
10x5mm

SI-1595
10x5mm

SI-1596
Bone Graft
2-6-6-8-10mm

SI-1597
10x5mm

SI-1598
10x5mm

SI-1599
10x

PERIOSTEALS ELEVATOR 10mm Hollow Handle



SI-1600
Raspator
BonexScraper
Length 165 mm
10 mm hollow
handle
double ended 6
mm/4 mm'



SI-1601
Raspator
Raspatory
Length 180 mm
10 mm hollow
handle
double ended 3
mm 3 mm'



SI-1602
Raspator
Raspatory
Length 180 mm
10 mm hollow
handle
double ended 3
mm 3 mm'



SI-1603
Raspator
Papilla
Length 185 mm
10 mm hollow
handle
double ended 3
mm 2 mm'



SI-1604
Raspator
Bone Scraper
Length 175 mm
10 mm hollow
handle
double ended



SI-1605
Raspator
Micro freer
Length 180 mm
6 mm solid grip
double ended



SI-1606
Raspator
Freer
Length 180 mm
6 mm solid handle
double ended
5mm 5 mm

PERIOSTEALS ELEVATOR 10mm Hollow Handle



SI-1607
Raspator
Freer
10 mm hollow
high quality
medical
Double ended

SI-1608
Raspator
Molt
Length 180 mm
8 mm solid grip
Double ended

SI-1609
Raspator
micro Freer
Length 185 mm
10 mm hollow
handle
Double ended

SI-1610
Raspator
micro Freer
Length 185 mm
10 mm hollow
handle
Double ended

SI-1611
Raspator
micro chisel
Length 170 mm
10 mm hollow
handle
Double ended

SI-1612
Raspator
micro chisel
Length 170 mm
10 mm hollow
handle
Double ended

SI-1613
Raspator
micro chisel
Length 175 mm
8 mm solid
handle
Double ended
3.5 mm 3.5mm

PERIOSTEALS ELEVATOR 10mm Hollow Handle



SI-1614
Raspator
micro Molt
Length 175 mm
corrosion resistant

SI-1615
Raspator Chipping
Chisel 24G and 7.5 mm
length about 175 mm
corrosion resistant

SI-1616
Raspator Molt
5.0 mm and 4.0 mm
length about 175 mm
corrosion resistant

SI-1617
Rasp Molt M9
Length 175 mm
corrosion resistant

SI-1618
Raspator Wautier
Length 175 mm
corrosion resistant

SI-1619
Raspator Benque
Length 170 mm
worker
5 mm 7mm
corrosion resistant

SI-1620
Raspator Nordberg
Length 180 mm
8 mm hollow handle
double ended
6 mm 4 mm

Hollow Handle 10mm Length 170mm



TISSUE MEMBRANE FORCEPS



AE-1630
Tissue Graft Forceps
Serrated Teeth



SI-1631
Self Locking
Grooved



SI-1632
Tissue Graft
forceps serrated
teeth



SI-1633
Titanium
Locking Plier



SI-1634
Titanium
Locking Plier



AE-1635
Adson Carbide
Tipped Tissue Forceps



AE-1636
Micro Adson 1x2
Tissue Forceps 12.1cm
1x2 Teeth 0.5mm tips



AE-1637
Adson 7x7
Tissue Forceps
7x7 Teeth 12cm



AE-1638
Adson 7x7 curved
Tissue Forceps
7x7 Teeth 12cm



AE-1639
Suture Pliers
Suture Forceps with
1.6mm Ring Tip



AE-1640
Debakey Atraumatic
Tissue Forceps Straight
Serrated Teeth
16cm



AE-1641
Debakey Atraumatic
Tissue Forceps Curved
Serrated Teeth
16cm

TISSUE MEMBRANE FORCEPS



AE-1642
Russian Long
Tissue Forceps
Serrated Tip
20 cm



AE-1643
Russian
Tissue Forceps
Serrated Tip
15 cm



SI-1644
Tissue Punch Forceps
to remove tissue
After Using Punch



SI-1645
13.5cm
narrow



SI-1646
10.5 cm



SI-1647
13.5 cm
wide



AE-1648
Pin Tweezers
with or without lock
length 15 cm
with lock



SI-1649
without lock



SI-1650
Haltepinzette with lock
15,5 cm for pins towards
the tips etc



SI-1651
Carry Forceps
15cm



SI-1652
Carry Forceps
15cm



SI-1653
Carry Forceps
15cm

OSTEOTOME CONVEXES / CONCAVE / CONIQUES

CONVEXES



AE-1654
SET OF 5 PIECES

CONCAVES




SI- 1655
SET OF 5 PIECES

CONIQUES



AE-1656
SET OF 5 PIECES



	CONVEXES STRAIGHT	1654	2,7 Ømm	3,2 Ømm	3,7 Ømm	4,2 Ømm	5 Ømm
	CONCAVES STRAIGHT	1655	2,7 Ømm	3,2 Ømm	3,7 Ømm	4,2 Ømm	5 Ømm
	CONIQUES STRAIGHT	1656	2,7 Ømm	3,2 Ømm	3,7 Ømm	4,2 Ømm	5 Ømm


OSTEOTOME CONVEXES / CONCAVE / CONIQUES




AE-1657
SET OF 5 PIECES

AE-1658
SET OF 5 PIECES

SI-1659
SET OF 5 PIECES



	CONVEXES CURVED	1657	2,7 Ømm	3,2 Ømm	3,7 Ømm	4,2 Ømm	5 Ømm
	CONCAVES CURVED	1658	2,7 Ømm	3,2 Ømm	3,7 Ømm	4,2 Ømm	5 Ømm
	CONIQUES CURVED	1659	2,7 Ømm	3,2 Ømm	3,7 Ømm	4,2 Ømm	5 Ømm



SI-1660
Int. 2mm
Ext. 3mm



SI-1661
Int. 3mm
Ext. 3.5mm



SI-1662
Int. 4mm
Ext. 5mm



SI-1663
Int. 5mm
Ext. 6mm



SI-1664
Int. 6mm
Ext. 7mm



SI-1665
Int. 7mm
Ext. 8mm



SI-1666
Int. 8mm
Ext. 9mm



SI-1667
Int. 9mm
Ext. 10mm



SI-1668
Int. 10mm
Ext. 11mm

4mm 5mm 6mm 7mm



SI-1669
Tissue Punch with pin
Set of four pieces



SI-1670
Tissue Punch
Set of Five Pieces



SI-1671
Tissue Punch
Set of Five Pieces



SI-1672
Trephine Drill
Set of Eight Pieces



2.0mm 2.5mm 2.8mm

SI-1673
Short Implant Drills

3.2mm 3.65mm 4.2mm 5.2mm

SI-1674
Long Implant Drills



SI-1675
Trephine Drill
Set of Eight Pieces


BONE SCRAPER STRAIGHT & CURVED



SI-1676
Bone Scraper
Straight

SI-1677
Bone Scraper
Curved

TISSUE PUNCH STRAIGHT & CURVED



SI-1678
Tissue Paper
Straight

SI-1685
Tissue Paper
Curved



SI-1692
Bone Mill



SI-1693
Basins for Bone



SI-1694
BONE CRUSHER

KEYS FOR IMPLANT SCREW



SI-1695



SI-1696



SI-1697



SI-1698

TORQUE WRENCH & RATCHET



AE-1699 Star Torque Wrench 6.35mm Universal



SI-1700 Round Ratchet Wrench 4.0mm Universal



SI-1701 Hex Driver With 4.0mm Driver



AE-1702 Star Torque Wrench 6.35mm Universal



AE-1703 Hex Torque Wrench 4.00mm Universal



AE-1704 Star Torque Wrench 6.35mm Universal



AE-1705 Star Torque Wrench 6.35mm Universal



AE-1706 Star Torque Wrench 6.35mm Universal



AE-1707 Star Torque Wrench 6.35mm Universal



AE-1708 Star Torque Wrench 6.35mm Universal



SI-1709 Star Ratchet Wrench 6.35mm Universal



SI-1710 Star Ratchet Wrench 6.35mm Universal



SI-1711 Star Ratchet Wrench 6.35mm Universal



SI-1712 Star Ratchet Wrench 4.0mm Universal



SI-1713 Round Ratchet Wrench 4.0mm Universal



SI-1714 Round Ratchet Wrench 3.35mm Universal



SI-1715 Star Ratchet Wrench 4.0mm Universal



SI-1716 Star Ratchet Wrench 6.35mm Universal



BONE EXPANDER SCREW & SURGICAL KITS



SI-1722

Bone Universal Kit
 Drill: 2mm, 2.8mm, 3.2mm, 3.65mm,
 4.2mm, 4.5mm, 4.8mm, 5.2mm
 Trephines Bit:
 3.0mm, 4.0mm, 5.0mm, 6.0mm,
 7.0mm, 8.0mm, 9.0mm, 10mm.
 Driver:
 2.5mm, 2.5mm Hand Pcs
 1.25mm, 1.25mm Hand Pcs
 Short Driver 2.5mm
 Long Driver 2.5mm
 1.25mm, Short Hex Driver Hand Pcs
 1.25mm, Long Hex Driver Hand Pcs
 Torque Wrench 4.0mm



SI-1723

Bone Expander Kit
 Expander: 2.6mm, 3.0mm,
 3.4mm, 3.8mm, 4.3mm
 Ratchet Driver
 Parallel Pin
 Pilot Drill
 Drill 1.4mm
 Drill Driver 4mm
 Hand Driver 4mm
 Ratchet 4.0mm



SI-1724

Bone Universal Kit
 Drill: 2mm, 2.8mm, 3.2mm, 3.65mm,
 4.2mm, 4.5mm, 4.8mm, 5.2mm
 Trephines Bit:
 3.0mm, 4.0mm, 5.0mm, 6.0mm,
 7.0mm, 8.0mm, 9.0mm, 10mm.

Bone Compression:
 2mm, 3.3mm, 2.8mm, 4.2mm, 3mm, 3.5mm
 Driver:
 2.5mm, 2.5mm Hand Pcs
 1.25mm, 1.25mm Hand Pcs
 Short Driver 2.5mm
 Long Driver 2.5mm
 1.25mm, Short Hex Driver Hand Pcs
 1.25mm, Long Hex Driver Hand Pcs
 Torque Wrench 6.35mm 1 Pcs



SI-1725

Conical Drills: 20 Pcs
 Tissue Punches:
 3.0mm, 3.50mm, 4.0mm, 4.5mm,
 5.0mm, 6.0mm, 6.50mm
 Bone Expander:
 1.6/3.8mm, 1.9/3.0mm, 2.1/3.6mm,
 2.6/3.8mm, 2.9/4.3mm
 Drill Adapter: 1 Pcs
 1.25mm, 2.50mm Hex Driver Hand Pcs

Pilot Drill: 1.6mm
 Drill Bit: 1.4mm, 1.6mm, 1.8mm, 2.0mm, 2.2mm
 Ratchet Driver: 4mm/6.35, 1.25/2.50
 Rimer: 3.2/4.7mm, 3.7/4.7mm
 Parallel Pin: 4.5mm, 2.5mm
 Double step Drill: 2.4/2.8mm
 Depth Gauge: 1 Pcs
 Torque Wrench 6.35mm Universal: 1 Pcs

BONE INJECTORS, COLLECTORS & SUCTION TUBE



SI-1726
2.5mm
3.5mm
170mm

SI-1727
4.5mm
5.5mm
170mm

SI-1728
3.5mm
4.5mm
170mm

SI-1729
6.0mm
7.0mm
170mm



SI-1730
18cm
Straight/ Bend
I.D 2.5mm
O.D.3.5

SI-1731
18cm
Straight/ Bend
I.D 3.5mm
O.D.4.5

SI-1732
18cm
Straight/ Bend
I.D 4.5mm
O.D.5.5

SI-1733
18cm
Straight/ Bend
I.D 6.0mm
O.D.7.0



SI-1734
17.5cm

SI-1735
18cm

SI-1736
coarse
5mm

SI-1737
Fine
9mm



SI-1738
Suction
tubes
180mm

SI-1739
Suction
tubes
180mm

SI-1740
Suction
tubes
180mm



SI-1742
3 mm
length 180mm


SI-1743
4 mm
length 180mm

SI-1744
5 mm
length 180mm

SI-1741
Fig. 1
Frazier
15cm -2 mm

PERFORATED WITH RIM LOCK

Set of 06 Pcs




SI-1745

Set of 6 PCS
Upper S.M.L
Lower S.M.L

PERFORATED WITH RIM LOCK

Set of 8 Pcs



SI-1746

Set of 8 PCS
Upper XS S.M.L
Lower XS S.M.L

SI-1747

Set of 8 PCS
Upper S.M.L,XL
Lower S.M.L,XL

SI-1748

Set of 10 PCS
Upper XS,S.M.L,XL
Lower XS,S.M.L,XL

Set of 6 Pcs



SI-1749
Set of 6 PCS
Upper S.M.L
Lower S.M.L

SOLID WITH RIM LOCK
Set of 8 Pcs




SI-1750
Set of 8 PCS
Upper XS S.M.L
Lower XS S.M.L

SI-1751
Set of 8 PCS
Upper S.M.L.XL
Lower S.M.L.XL

SI-1752
Set of 10 PCS
Upper XS,S.M.L,XL
Lower XS,S.M.L,XL

PERFORATED FULL DENTURE

Set of 8 Pcs



SI-1753

Set of 8 PCS
Upper S.M.L.XL
Lower S.M.L.XL

SI-1754

Set of 8 PCS
Upper S.M.L
Lower S.M.L

SOLID FULL DENTURE

Set of 8 Pcs




SI-1755

Set of 8 PCS
Upper S.M.L.XL
Lower S.M.L.XL


SI-1756

Set of 8 PCS
Upper S.M.L
Lower S.M.L

Perforated Without Rim Lock



Solid With Rim Lock






SI-1793
Instruments Cassettes
5 Pieces, 7 Pieces, 10 Pieces & 20 Pieces



SI-1794
Instruments Cassettes
5 Pieces, 7 Pieces, 10 Pieces & 20 Pieces



SI-1795
Instruments Cassettes
7 Pieces & 10 Pieces



SI-1796
Gracey Curette Set Cassettes



SI-1797
182 x 150 x 20 mm
10 Instruments



SI-1798
180 x 85 x 15 mm
5 Instruments



SI-1799
288 x 182 x 35 mm
20 Instruments



SI-1800
288 x 182 x 35 mm
20 Instruments




SI-1801
182 x 145 x 35 mm
10 Instruments




SI-1802
182 x 72 x 35 mm
5 Instruments



SI-1803
182 x 145 x 35 mm
9 Instruments



SI-1804
Gracey Curette Set Stand



SI-1805
Sinus Lift Set Stand



KIDNEY TRAY WITH COVER
(English Pattern)

SI-1806 150mm
SI-1807 200mm
SI-1808 240mm
SI-1809 300mm



KIDNEY TRAYS

With lid available

With lid available

SI-1810 6"
SI-1811 8"
SI-1812 10"
SI-1813 12"

SI-1814 6"
SI-1815 8"
SI-1816 10"
SI-1817 12"



Deep

SI-1818 2.25"- 2"
SI-1819 4.75"- 2"
SI-1820 5.5"- 2.5"
SI-1821 6.25"- 2.75"
SI-1822 7" 3.5"



Sterilizing Drum

Art #	mm	mm
SI-1823	80	80
SI-1824	120	120
SI-1825	160	120
SI-1826	190	140
SI-1827	190	200
SI-1828	240	120
SI-1829	240	160
SI-1830	240	240
SI-1831	290	120
SI-1832	290	160
SI-1833	290	240
SI-1834	340	160
SI-1835	340	240
SI-1836	390	145
SI-1837	390	200
SI-1838	390	290



SI-1839
Cotton Holder



SI-1840
Spirit Lamp



Endo Box (72 Holes)

SI-1841
Endodontic Box



Forcep Jar

SI-1842 25 x 80mm
SI-1843 30 x 100mm
SI-1844 100 x 200mm
SI-1845 75 x 173mm



SI-1846 16" x 8"
SI-1847 8" x 10"
SI-1848 10" x 12"
Without Lid



SI-1849
Instruments Box



SI-1850
Scaler Tray
100 x 200 mm

Surgimpex & Co.



E-MAIL

Info@surgimpex.com



PHONE

+923316163388



WEBSITE

www.surgimpex.com



ADDRESS

Defence Road Mustafa Town
S I E Sialkot Pakistan

